



Medications:

Epinephrine:
 0.01 mg/kg
 q 3-5 minutes
 (every other cycle)

Amiodarone
 5mg/kg
 (up to 3 doses)

or

Lidocaine
 1mg/kg loading dose
 (follow with drip)

ETCO2:

Normal 35-45 mmHg
 CPR: >10mmHg
 Spike = ROSC?

@1:45:

1. Pre-charge Defibrillator
2. Check pulse w/ CPR
3. Prepare to switch compressors

@2:00:

1. Stop CPR/Analyze/Shock?
2. Resume CPR
3. Push Med
4. Prep Next Med (closed loop)

ROSC

V - ITALS
 (O2: 94-99%, BP>5%)

A - IRWAY

R - ESPONSIVENESS

M - ONITOR (EKG)

I - V x2, LABS

N - EURO
 (CT, MRI, EEG)

T - TM / TRANSFER
 (5 days normothermia
 or 2/3 @32-34 C)
 (OR/ICU)